

---

# **Time Tagger Documentation**

***Release 1.0.4***

**Swabian Instruments**

**Nov 25, 2016**



## CONTENTS

<b>1</b>	<b>Installation instructions</b>	<b>1</b>
1.1	Windows users . . . . .	1
1.2	Linux users . . . . .	2
<b>2</b>	<b>Getting Started</b>	<b>3</b>
2.1	Windows users . . . . .	3
<b>3</b>	<b>Hardware</b>	<b>5</b>
3.1	Input channels . . . . .	5
3.2	Data connection . . . . .	5
3.3	Status LEDs . . . . .	5
3.4	Test signal . . . . .	6
3.5	Virtual channels . . . . .	6
3.6	Synthetic input delay . . . . .	6
3.7	Synthetic dead time . . . . .	6
3.8	Event filter . . . . .	6
3.9	Bin equilibration . . . . .	7
3.10	Overflows . . . . .	7
3.11	General purpose IO (available upon request) . . . . .	7
<b>4</b>	<b>Software Overview</b>	<b>9</b>
4.1	Web application and JSON-RPC interface . . . . .	9
4.2	Precompiled libraries and high level language bindings . . . . .	9
4.3	C++ API . . . . .	9
<b>5</b>	<b>Application Programmer's Interface</b>	<b>11</b>
5.1	Overview . . . . .	11
5.2	Organization . . . . .	12
5.3	The TimeTagger class . . . . .	12
5.4	Measurement Classes . . . . .	14
<b>6</b>	<b>Revision History</b>	<b>21</b>
6.1	V1.0.4 - 24.11.2016 . . . . .	21
6.2	V1.0.2 - 28.07.2016 . . . . .	21
6.3	V1.0.0 . . . . .	22



## INSTALLATION INSTRUCTIONS

### 1.1 Windows users

#### 1.1.1 Requirements

##### Operating System

The installer provided supports Windows 7 and Windows 10 (32 and 64 bit).

##### Python

If you want to run the Python **quickstart guide**, you need the following additional third-party software packages (available for free).

- **Python 2.7.x** or **Python 3.5.x** installed along with **ipython**
- the following Python modules: **numpy**, **matplotlib (pylab)**, **time**, **cPickle**

On Windows it is easiest to install a **Python distribution** that contains all required packages by default. The **anaconda** or **canopy** Python distributions will cover this. If you don't have Python installed, please download and install now either of these Python distributions. Note that **anaconda** is available as **Python 3.5.x**.

---

**Note:** If you are using a fresh install of the **canopy** distribution, please make sure to run the graphical interface once from the start menu to have your Python paths set up correctly.

---

##### Installation

---

**Note:** If you intend to use the *TimeTagger 20* with the Python programming language, please make sure that you install Python first as described above. This ensures that the installer detects your Python installation automatically.

---

1. Download and run the most recent *TimeTagger 20* windows installer.

Connect the *TimeTagger 20* to your computer with the USB cable.

**Caution:** Wait until windows has recognized the device and finished installing the driver.

You should now be ready to use your *TimeTagger 20*.

## 1.2 Linux users

The TimeTagger can be used under Linux with C++ or Python. There is no installer package yet, so please contact us if you want to use the TimeTagger under Linux.

## GETTING STARTED

### 2.1 Windows users

1. Make sure the *TimeTagger 20* software and a Python distribution are installed and the *TimeTagger 20* connected to your computer (see previous section).
2. Open a command shell and `cd` to the `.\examples\python\quickstart` folder under your *TimeTagger 20* installation directory
3. Start an **ipython** shell with plotting support by entering `ipython --pylab`
4. Run the **quickstart.py** script by entering `run quickstart`

The script demonstrates a selection of the features provided by the *TimeTagger 20* programming interface and runs some example measurements using the built in test signal generator and plots the result.

You are encouraged to open and read the `quickstart.py` file in an editor to see what it is doing.

Among others, the script will...

1. Create an instance called 'tagger' that represents the device.
2. Start the built in test signal (~0.8 MHz square wave) and apply it to channels 0 and 1
3. Create a time trace of the click rate on the first two channels, let it run for a while and plot the result.
4. Create coarse and fine cross correlation measurements. The coarse measurement shows characteristic peaks at integer multiples of the inverse frequency of the test signal. The fine measurement demonstrates the < 60 ps time resolution.
5. Show you how to create virtual channels, use synchronization, event filter and control the input trigger level.

#### 2.1.1 Where to go from here

To learn more about the *TimeTagger 20* you are encouraged to consult the following resources.

1. If you have not done so already, have a look at the Python script you just run.
2. More details about the software interface are covered by the API documentation in the subsequent section
3. You are also encouraged to study the C++ source code provided in the *TimeTagger 20* install directory
4. Code examples in `.\examples\python\traits` show you how to quickly implement an interactive graphical interface in python

---

**Note:** The python GUI examples, additionally require the following python packages installed on your system: *traits*, *traitsui*, *chaco*, *pyface enable*.

---

### 2.1.2 DLL/wrappers for .net (e.g. Matlab, LabView)

We provide for the TimeTagger a .net class Library (32 and 64 bit) which can be used to access the TimeTagger from high level languages `.\dotNet\TTCSharpxx.dll`.

The following is important to note:

- Namespace: `SwabianInstruments.TimeTagger`
- static functions (e.g. `create an instance of a TimeTagger`) are accessible via `SwabianInstruments.TimeTagger.TT`

### 2.1.3 LabView

A simple correlation measurement is provided in `.\examples\LabView\` for LabView 2014. The requirements for Using .NET assemblies in LabVIEW can be found here: \* Link <[http://zone.ni.com/reference/en-XX/help/371361L-01/lvconcepts/net\\_defaults/](http://zone.ni.com/reference/en-XX/help/371361L-01/lvconcepts/net_defaults/)>

---

**Note:** LabView most likely is installed as a 32 bit version so include the 32 bit .net library to your project

---

### 2.1.4 Matlab

Wrapper classes are provided for Matlab so that native Matlab variables can be used.

The TimeTagger toolbox is automatically installed during the setup. If TimeTagger is not available in your Matlab environment try to reinstall the toolbox from `.\MatlabWrapper\TimeTagger Matlab.mltbx`.

The following changes in respect to the .net library have been made:

- static functions are available through the `TimeTagger` class
- all classes except for the `TimeTagger` class itself have a `TT` prefix (e.g. `TTCounter`) to not collide with any variables/classes in your Matlab environment

An example how to use the TimeTagger with Matlab can be found in `.\examples\MatLab\`.



## HARDWARE

### 3.1 Input channels

The *TimeTagger 20* has 8 SMA connectorized input channels numbered 0 to 7 throughout this document. The electrical characteristics are tabulated below. Both rising and falling edges are detected on the input channels. On the software level, rising edges correspond to channel numbers 0 to 7 and falling edges correspond to respective channel numbers 8 to 15. Thereby, you can treat rising and falling edges in a fully equivalent fashion.

#### 3.1.1 Electrical characteristics

Property	Value
Termination	50 $\Omega$
Input voltage range	0 to 5 V
Trigger level range	0 to 3.3 V
Minimum signal level	~50 mV
Minimum pulse width	~1 ns

### 3.2 Data connection

A USB connection is used for data and power supply. Please ensure that the USB port is capable of providing the full specified current (500 mA). A USB 2.0 data connection is required for reasonable performance. Operating the device via a USB hub is strongly discouraged. The *TimeTagger 20* can stream about 5 M tags per second.

### 3.3 Status LEDs

The *TimeTagger 20* has two LEDs showing status information. A green LED turns on when the USB power is connected. An RGB LED shows the information tabulated below.

green	firmware loaded
blinking green-orange	time tags are streaming
red flash (0.1 s)	an overflow occurred
continuous red	repeated overflows

## 3.4 Test signal

The *TimeTagger 20* has a built in test signal generator that generates a square wave with a frequency in the range 0.9 to 1.0 MHz. You can apply the test signal to any input channel instead of the external input. This is useful for testing, calibrating and setting up,

## 3.5 Virtual channels

The architecture allows you to create virtual channels, e.g., you can create a new channel that represents the sum of two channels (logical OR), or coincidence clicks of two channels (logical AND).

## 3.6 Synthetic input delay

You can introduce for each channel an input delay. This is useful e.g. to compensate for propagation delay in cables of unequal length, if the relative timing between two channels is important. The input delay can be set individually for rising and for falling edges.

## 3.7 Synthetic dead time

You can introduce for each channel a synthetic dead time. This is useful when you want to suppress consecutive clicks that are closely separated, e.g., to suppress after-pulsing of avalanche photo diodes or to suppress too high data rates. The dead time can be set individually for rising and for falling edges.

## 3.8 Event filter

The event filter allows you to decrease the time tag rate without losing those time tags that are relevant to your application. In a typical use case you have a high frequency signal applied to at least one channel. Examples include fluorescence lifetime measurements or optical quantum information and cryptography where you want to capture synchronization clicks from a high repetition rate excitation laser.

To reduce the data rate you discard all synchronization clicks, except those that succeed one of your low rate detector clicks, thereby forming a reduced time tag stream. The reduced time tag stream is processed by the software in the exact same fashion as the full time tag stream.

This feature is enabled by the event filter. All channels on your Time Tagger 20 are fully equivalent. You can specify which channels are filtered and which channels are used as triggers that enable transmission of a subsequent tag on the filtered channels.

**Caution:** The time resolution of the event filter is 6 ns. To ensure deterministic filter logic, the physical time difference between the filtered channels and triggered channels must be larger than 9 ns. Note also that software defined input delays as set by the method `setInputDelay()` do not apply to the event filter logic.

## 3.9 Bin equilibration

Discretization of electrical signals is never perfect. In time-to-digital conversion, this is manifest as small differences (few ps) of the bin sizes inside the converter that even varies from chip to chip. This imperfection is inherent to any time-to-digital conversion hardware. It is usually not apparent to the user. However, when correlations between two channels are measured on short time scales you might see this as a weak periodic ripple on top of your signal. We balance this on the software level at the cost of a decrease of the time resolution by  $\sqrt{2}$ . This feature is enabled by default. If your application requires time resolution down to the jitter limit, you can disable this feature.

## 3.10 Overflows

The *TimeTagger 20* is capable of streaming on average about 5 million tags per second. Higher data rates for short times will be buffered internally so that no overflow occurs. This internal buffer is limited so that on continuously higher data rates, data loss occurs and parts of the time tags are lost. The hardware allows you to check with `timeTagger.getOverflows()` whether an overflow condition has occurred. If no overflow is returned, you can be sure that every time tag is received.

---

**Note:** When overflows occur, Time Tagger 20 will still produce valid chunks of data and discard the invalid tags inbetween. Your measurements data may still be valid, albeit, your acquisition time likely goes up.

---

## 3.11 General purpose IO (available upon request)

The device is ready to be equipped with up to four SMA connectorized general purpose IO ports and an external clock input or output. These can be used to implement custom features such as special fast input or output triggers, enable / disable gates, software controllable input and output lines, etc.. Please contact us for custom designs.



## SOFTWARE OVERVIEW

The heart of the *TimeTagger 20* software is a multi-threaded driver that receives the time tag stream and feeds it to all running measurements. Measurements are small threads that analyze the time tag stream each in their own way. For example, a count rate measurement will extract all time tags of a specific channel and calculate the average number of tags received per second, a cross-correlation measurement will compute the cross-correlation between two channels, typically by sorting the time tags in histograms, etc.. This is a powerful architecture that allows you to perform any thinkable digital time domain measurement in real time. You have several choices to use this architecture.

### 4.1 Web application and JSON-RPC interface

---

**Note:** This feature is currently under development. Please contact us if you would like to use the TimeTagger in this way.

---

The easiest way of using the *TimeTagger 20* is via a web application that allows you to interact with the hardware from a web browser on your computer or a tablet. You can create measurements and get live plots, save and load the acquired data from within a web browser. In addition, you can also access and remote control the web application via a JSON-RPC interface.

### 4.2 Precompiled libraries and high level language bindings

We have implemented a set of typical measurements including count rates, auto correlation, cross correlation, fluorescence lifetime imaging (FLIM), etc.. For most users, these measurements will cover all needs. These measurements are included in the C++ API and provided as precompiled library files. To make using the TimeTagger even easier, we have equipped these libraries with bindings to higher level languages (Python, Matlab, LabView, .net in general) so that you can directly use the TimeTagger from this language. With these APIs you can easily start a complex measurement from a higher level language with only two lines of code. To use one of these APIs, you have to write code in the high level language of your choice. Refer to chapter [Getting Started](#) and [Application Programmer's Interface](#) if you plan to use the TimeTagger in this way.

### 4.3 C++ API

The underlying software architecture is provided by a C++ API that implements two classes: one class that represent the TimeTagger and one class that represents a base measurement. On top of that, the C++ API also provides all predefined measurements that are made available by the web application and high level language bindings. To use this API, you have to write and compile a C++ program.



## APPLICATION PROGRAMMER'S INTERFACE

### 5.1 Overview

The API provides methods to control the hardware and to create measurements that are hooked onto the time tag stream. It is written in C++ but wrapper classes for higher level languages (python, Matlab, LabView) are provided, such that the C++ API can directly be used in your application, in a way that is equivalent to the C++ classes. The API includes a set of typical measurements that will most likely cover your needs. Implementation of custom measurements is based on subclassing from a C++ base class and thus is only available in the C++ API.

#### 5.1.1 API documentation

The API documentation in this manual gives a general overview how to use the *TimeTagger 20*.

#### 5.1.2 Examples

Often the fastest way to learn how to use an API is by means of examples. Please see the `\examples` subfolder of your *TimeTagger 20* installation for examples.

#### 5.1.3 Units

Time is measured in ps since device startup and represented by 64 bit integers. Note that this implies that the time variable will rollover once after about 0.83 years. This will most likely not be relevant to you unless you plan to run your software continuously over one year and are taking data at the instance when the rollover is happening.

#### 5.1.4 Channel Numbers

You can use the *TimeTagger 20* to detect both rising and falling edges. Throughout the software API, the rising edges are represented by channels 0 to 7 and the falling edges are represented by channel numbers 8 to 15. Virtual channels will obtain numbers from 16 onwards.

#### 5.1.5 Undefined Channels

There might be the need to leave a parameter undefined when calling a class constructor. Depending on the programming language you are using you pass an undefined channel via the static constant `CHANNEL_UNDEFINED` which can be found in the `TT` class for .net and in the `TimeTagger` class in Matlab.

## 5.2 Organization

The API contains a *small* number of **classes** which you instantiate in your code. These **classes** are summarized below.

### 5.2.1 Hardware

**TimeTagger** Represents the hardware and provides methods to control the trigger levels, input delay, dead time, event filter and test signals.

### 5.2.2 Virtual Channels

**Combiner** Combines two channels into one

**Coincidence** Detects coincidence clicks on two or more channels within a given window

**Coincidences** Detects coincidence clicks on multiple channel groups within a given window

### 5.2.3 Measurements

**Correlation** auto- and cross-correlation.

**CountBetweenMarkers** Counts tags on one channel where the bins are determined by triggers on one or two other channels. Uses static buffer output. Use this to implement a gated counter, a counter synchronized to an external sampling clock, etc.

**Counter** Counts clicks on one or more channels with a fixed binwidth and circular buffer output.

**Countrate** Average tag rate on one or more channels.

**FLIM** Fluorescence lifetime imaging.

**Iterator** Base class for implementing custom measurements.

**Histogram** A simple histogram of time differences. This can be used e.g. to measure lifetime.

**Scope (experimental)** Detects the rising and falling edges on a channel to visualize the incoming signals similar to an ultrafast digitizer

**StartStop** Accumulates a histogram of time difference between pairs of tags on two channels. Only the first stop tag after a start tag is considered. Subsequent stop tags are discarded. The Histogram length is unlimited.

**TimeDifferences** Accumulates the time differences between tags on two channels in one or more histograms. The sweeping through histograms is optionally controlled by one or two additional triggers.

## 5.3 The TimeTagger class

This class provides access to the hardware and exposes methods to control hardware settings. Behind the scenes it opens the USB connection, initializes the device and receives the time tag stream. Every measurement requires an instance of the TimeTagger class to which it will be associated. In a typical application you will perform the following steps:

1. create an instance of TimeTagger
2. use methods on the instance of TimeTagger to adjust the trigger levels
3. create an instance of a measurement passing the instance of TimeTagger to the constructor



You can use multiple TimeTaggers on one computer simultaneously. In this case, you usually want to associate your instance of the TimeTagger class to a physical TimeTagger. To implement this in a bullet proof way, TimeTagger instances must be created with a factory function called 'createTimeTagger'. The factory function accepts the serial number of a physical TimeTagger as a string argument (every TimeTagger has a unique hardware serial number). The serial number is the only argument that can be passed. If an empty string or no argument is passed, the first detected TimeTagger will be used. To find out the hardware serial number, you can connect a single TimeTagger, open it and use the 'getSerial' function described below.

The TimeTagger class contains a small number of methods to control the hardware settings that are summarized below.

### 5.3.1 Methods

**reset()** reset the TimeTagger to the startup state

**setTriggerLevel(unsigned int channel, double voltage)** set the trigger level of an input channel in volts

**double getTriggerLevel(unsigned int channel)** return the trigger level of an input channel in volts

**setInputDelay(unsigned int channel, long long delay)** set the input delay of a channel in picoseconds

**long long getInputDelay(unsigned int channel)** return the input delay of a channel in picoseconds

**setFilter(bool state)** Deprecated since version V1.0.3: use `setConditionalFilter()` instead.

enable or disable the event filter on the FPGA board. If the filter is active, tags on channel 7 are only transmitted if they were immediately preceded by a tag on channel 0.

**bool getFilter()** Deprecated since version V1.0.3: use `getConditionalFilter()` instead.

returns true if the event filter on the FPGA board is enabled

**setConditionalFilter(<vector int> trigger, <vector int> filter)** New in version V1.0.3: activates or deactivates the event filter.

Time tags on the filter channels are discarded unless they were preceded by a time tag on one of the trigger channels.

**setNormalization(bool state)** New in version V1.0.3: enables or disable Gaussian normalization of the detection jitter. Enabled by default.

**bool getNormalization()** returns true if Gaussian normalization is enabled.

**long long setDeadTime(unsigned int channel, long long deadtime)** sets the dead time of a channel in picoseconds. The requested time will be rounded to the nearest multiple of the clock time. The deadtime will also be clamped to device specific limitations. As the actual deadtime will be altered, the real value will be returned.

**long long getDeadTime(unsigned int channel)** returns the dead time of a channel in picoseconds

**setTestSignal(unsigned int, bool state)** This will connect or disconnect the channel with the on-chip uncorrelated signal generator.

**bool getTestSignal(unsigned int channel)** returns true if the internal test signal is activated on the specified channel

**string getSerial()** returns the hardware serial number

**long long getOverflows()** returns the number of overflows that occurred since startup

**long long getOverflowsAndClear()** returns the number of overflows that occurred since startup and sets them to zero

**clearOverflows()** set the overflow counter to zero

**sync()** ensure that all hardware settings such as trigger levels, channel registrations, etc., have propagated to the FPGA and are physically active and synchronize the TimeTagger pipeline, so that all tags arriving after a sync call were actually produced after the sync call. The sync function waits until all historic tags in the pipeline are processed.

### 5.3.2 Debug Methods

**int getBoardModel()** returns the hardware type

**registerChannel(unsigned int channel)** enable transmissions of time tags on the specified channel

**unregisterChannel(unsigned int channel)** disable transmissions of time tags on the specified channel

**unsigned int getChannels()** returns the number of registered channels

**autoCalibration(bool verbose=true)** run an auto calibration of the tagger hardware using the built in test signal

**2D array long long getDistributionCount()** returns the calibration data represented in counts

**2D array long long getDistributionPSec()** returns the calibration data in picoseconds

**long long getPsPerClock()** returns the the duration of a clock cycle in picoseconds

## 5.4 Measurement Classes

The library includes a number of common measurements that will be described in this section. All measurements are derived from a base class called 'Iterator' that is described further down. As the name suggests, it uses the *iterator* programming concept.

All measurements provide a small number of methods to start and stop the execution and to access the accumulated data. The methods are summarized below.

### 5.4.1 Methods common to all Measurements

**getData()** Returns the data accumulated up to now. The returned data can be a scalar, vector or array, depending on the measurement.

**clear()** reset the accumulated data to an array filled with zeros

**start()** start data acquisition

**stop()** stop data acquisition

**Attention:** All measurements start accumulating data immediately after their creation.

In a typical application you will perform the following steps:

1. create an instance of a measurement, e.g. ~a countrate on channel 0
2. wait for some time
3. retrieve the data accumulated by the measurement up to now by calling the 'getData' method.

The specific measurements are described below.

### 5.4.2 Correlation

Accumulates time differences between clicks on two channels into a histogram, where all ticks are considered both as start and stop clicks and both positive and negative time differences are considered.

## Arguments

**tagger** <reference> reference to a time tagger

**channel 1** <int> first channel

**channel 2** <int> second channel

**binwidth** <longlong> binwidth in ps

**n\_bins** <int> the number of bins in the resulting histogram

## Methods

**getData()** returns a one-dimensional array of size  $2*n\_bins+1$  containing the histograms.

**getIndex()** returns a vector of size 'n\_bins' containing the time bins in ps.

**clear()** resets the array to zero.

**setMaxCounts()** set the maximum number of start tags accepted

**getCounts()** returns the number of start tags

**ready()** returns 'true' when the required number of start tags set by 'setMaxCounts' has been reached

### 5.4.3 CountBetweenMarkers

Count on a single channel. The bin edges between which counts are accumulated are determined by one or more hardware triggers. Specifically, the measurement records data into a vector of length 'n\_values' (initially filled with zeros). It waits for tags on the 'begin\_channel'. When a tag is detected on the 'begin\_channel' it starts counting tags on the 'click\_channel'. When the next tag is detected on the 'begin\_channel' it stores the current counter value as next entry in the data vector, resets the counter to zero and starts accumulating counts again. If an 'end\_channel' is specified, the measurement stores the current counter value and resets the counter when a tag is detected on the 'end\_channel' rather than the 'begin\_channel'. You can use this e.g., to accumulate counts within a gate by using rising edges on one channel as the 'begin\_channel' and falling edges on the same channel as the 'end\_channel'. The measurement stops when all entries in the data vector are filled.

## Arguments

**tagger** <reference> reference to a time tagger

**begin\_channel** <int> channel that triggers beginning of counting and stepping to the next value

**end\_channel** <int> channel that triggers end of counting

**n\_values** <int> number of values

## Methods

**getData()** returns an array of size 'n\_values' containing the acquired counter values.

**getIndex()** returns a vector of size 'n\_values' containing the time bins in ps.

**clear()** resets the array to zero and restarts the measurement.

**ready()** returns 'true' when the entire array is filled.

### 5.4.4 Counter

Time trace of the countrate on one or more channels. Specifically this measurement repeatedly counts tags on one or more channels within a time interval 'binwidth' and stores the results in a two dimensional array of size 'number of channels' times 'n\_values'. The array is treated as a circular buffer that is, all values in the array are shifted by one position when a new value is generated. The last entry in the array is always the most recent value.

#### Arguments

**tagger** <reference> reference to a time tagger

**channels** <vector int> channels used for counting tags

**binwidth** <longlong> binwidth in ps

**n\_values** <int> number of values

#### Methods

**getData()** returns an array of size 'number of channels' times 'n\_values' containing the current values of the circular buffer (counts in each bin).

**getIndex()** returns a vector of size 'n\_values' containing the time bins in ps.

**clear()** resets the array to zero and restarts the measurement.

### 5.4.5 Countrate

Measures the average countrate on one or more channels. Specifically, it counts tags on the specified channels and determines the time between the first tag since instantiation and the latest tag. The ratio of the number of tags and the time corresponds to the average countrate since the first tag.

#### Arguments

**tagger** <reference> reference to a time tagger

**channels** <vector int> channels used for counting tags

#### Methods

**getData()** returns the average countrate in counts per second.

**clear()** resets the accumulated counts to zero and uses the next incoming tag as the first tag.

### 5.4.6 FLIM

Fluorescence-lifetime imaging microscopy or FLIM is an imaging technique for producing an image based on the differences in the exponential decay rate of the fluorescence from a sample.

Fluorescence lifetimes can be determined in the time domain by using a pulsed source. When a population of fluorophores is excited by an ultrashort or delta pulse of light, the time-resolved fluorescence will decay exponentially.

This measurement implements a line scan in a FLIM image that consists of a sequence of pixels. This could either represent a single line of the image, or - if the image is represented as a single meandering line - this could represent the entire image.

This measurement is a special case of the more general 'TimeDifferences' measurement.

The measurement successively acquires n histograms (one for each pixel in the line scan), where each histogram is determined by the number of bins and the binwidth.

### Arguments

**tagger** <reference> reference to a time tagger

**click channel** <int> channel on which clicks are received

**start channel** <int> channel on which start clicks are received

**next\_channel** <int> channel on which pixel triggers are received

**binwidth** <longlong> binwidth in ps

**n\_bins** <int> number of bins in each histogram

**n\_pixels** <int> number of pixels

### Methods

**getData()** returns a two-dimensional array of size 'n\_bins' times 'n\_pixels' containing the histograms.

**getIndex()** returns a vector of size 'n\_bins' containing the time bins in ps.

**clear()** resets the array to zero.

**setMaxCounts()** set the maximum number of start tags accepted

**getCounts()** returns the number of start tags

**ready()** returns 'true' when the required number of start tags set by 'setMaxCounts' has been reached

## 5.4.7 Histogram

Accumulate time differences into a histogram. This is a simple multiple start, multiple stop measurement. This is a special case of the more general 'TimeDifferences' measurement. Specifically, the measurement waits for clicks on the 'start channel', and for each start click, it measures the time difference between the start click and all subsequent clicks on the 'click channel' and stores them in a histogram. The histogram range and resolution is specified by the number of bins and the binwidth specified in ps. Clicks that fall outside the histogram range are ignored. Data accumulation is performed independently for all start clicks. This type of measurement is frequently referred to as 'multiple start, multiple stop' measurement and corresponds to a full auto- or cross-correlation measurement.

### Arguments

**tagger** <reference> reference to a time tagger

**click channel** <int> channel on which clicks are received

**start channel** <int> channel on which start clicks are received

**binwidth** <longlong> binwidth in ps

**n\_bins** <int> the number of bins in the histogram

### Methods

**getData()** returns a one-dimensional array of size **n\_bins** containing the histogram.

**getIndex()** returns a vector of size 'n\_bins' containing the time bins in ps.

**clear()** resets the array to zero.

**setMaxCounts()** set the maximum number of start tags accepted

**getCounts()** returns the number of start tags

**ready()** returns 'true' when the required number of start tags set by 'setMaxCounts' has been reached

### 5.4.8 StartStop

A simple start-stop measurement. This class performs a start-stop measurement between two channels and stores the time differences in a histogram. The histogram resolution is specified beforehand (binwidth) but the histogram range is unlimited. It is adapted to the largest time difference that was detected. Thus all pairs of subsequent clicks are registered.

### Arguments

**tagger** <reference> reference to a time tagger

**click\_channel** <int> channel on which stop clicks are received

**start\_channel** <int> channel on which start clicks are received

**binwidth** <longlong> binwidth in ps

### Methods

**getData()** returns a one dimensional array containing the histogram

**getIndex()** returns a vector of the time bins in ps.

**clear()** resets the array to zero and restarts the measurement.

**ready()** returns 'true' when the entire array is filled.

### 5.4.9 TimeDifferences

A multidimensional histogram measurement optionally with up to three additional channels that control how to step through the indices of the histogram array. This is a very powerful and generic measurement. You can use it to record cross-correlation, lifetime measurements, fluorescence lifetime imaging and many more measurements based on pulsed excitation. Specifically, the measurement waits for a tag on the 'start\_channel', then measures the time difference between the start tag and all subsequent tags on the 'click\_channel' and stores them in a histogram. If no 'start\_channel' is specified, the 'click\_channel' is used as 'start\_channel' corresponding to an auto-correlation measurement. The histogram has a number of 'n\_bins' bins of binwidth 'binwidth'. Clicks that fall outside the histogram range are discarded. Data accumulation is performed independently for all start tags. This type of measurement is frequently referred to as 'single start, multiple stop' measurement and corresponds to a full auto- or cross-correlation measurement.

The data obtained from subsequent start tags can be accumulated into the same histogram (one-dimensional measurement) or into different histograms (two-dimensional measurement). In this way you can perform more general two-dimensional time-difference measurements. The parameter 'n\_histograms' specifies the number of histograms. After each tag on the 'next\_channel', the histogram index is incremented by one (and reset to zero after reaching the last valid index). You can also provide a synchronization trigger that resets the histogram index by specifying a 'sync\_channel'.

Typically, you will run the measurement indefinitely until stopped by the user. However, it is also possible to specify the maximum number of rollovers of the histogram index. In this case the measurement stops when the number of rollovers has reached the specified value. This means that both for a one-dimensional and for a two-dimensional measurement, it will measure until every histogram has seen the specified number of start tags.

## Arguments

**tagger** <reference> reference to a time tagger

**click channel** <int> channel that increments the count in a bin

**start channel** <int> channel that sets start times relative to which s on the click channel are measured

**next channel** <int> channel that increments the histogram index

**sync channel** <int> channel that resets the histogram index to zero

**binwidth** <longlong> binwidth in ps

**n\_bins** <int> number of bins in each histogram

**n\_histograms** <int> number of histograms

## Methods

**getData()** returns a two-dimensional array of size 'n\_bins' times 'n\_histograms' containing the histograms.

**getIndex()** returns a vector of size 'n\_bins' containing the time bins in ps.

**clear()** resets the array to zero.

**setMaxCounts()** set the maximum number of start clicks accepted

**getCounts()** returns the number of start clicks

**ready()** returns 'true' when the required number of start clicks set by 'setMaxCounts' has been reached

### 5.4.10 Dump

Dump the time tag stream to a file in a binary format.

## Arguments

**<str> filename** name of the file to dump to

## Methods

**stop()** stop the measurement





## REVISION HISTORY

### 6.1 V1.0.4 - 24.11.2016

#### Hardware changes

- extended event filter to multiple condition and filter channels
- improved jitter for channel 0
- channel delays might be different from previous version ( $< 1$  ns)

#### API changes

- new function setConditionalFilter allows for multiple filter and event channels (replaces setFilter)
- Scope class implements functionality to use the time tagger as a 50 GHz digitizer
- Coincidences class now can handle multiple coincidence groups which is much faster than multiple instances of Coincidence
- added examples for C++ and .net

Software changes \* improved GUI (Web application)

Bug fixes \* Matlab/LabVIEW does not require to have the Visual Studio Redistributable package installed

### 6.2 V1.0.2 - 28.07.2016

#### Major changes:

- LabView support including various example VIs
- Matlab support including various example scripts
- .net assembly / class library provided (32 and 64 bit)
- WebApp graphical user interface to get started without writing a single line of code
- Improved performance (multicore CPUs are supported)

#### API changes:

- reset() function added to reset a TimeTagger device to the startup state
- getOverflowsAndClear() and clearOverflows() introduced to be able to reset the overflow counter
- support for python 3.5 (32 and 64 bit) instead of 3.4

## 6.3 V1.0.0

initial release supporting python